

Alex Jacob

CHMG 146 Section 12

Date of Analysis 4/2/2021, 4/9/2021, 4/16/2021

### ***Lab 5: Alkalinity***

#### ***Abstract***

Alkalinity is often described as the capacity of water to resist acidification in terms of concentration. This concentration corresponds to the amount of acid which is added as a titrant. In this case, the milligrams per liter calcium concentrate is used. To determine this, three separate experiments are used. The first focuses on standardizing a sulfuric acid solution via titration with carbonate to find the average molarity of the sulfuric acid. The second focuses on determining the alkalinities of an unknown substance using titration and visual cues. The third focuses on determining the alkalinities of the same solution but using pH probes to measure the rate of change in pH then graphing to find the endpoints. The unknown sample was determined to be a solution of  $\text{CO}_3^{2-}$  and  $\text{HCO}_3^-$ . The total alkalinity of this solution was calculated to be 2833.33 mg of  $\text{CaCO}_3$  / L using Part 2's method and 3125 mg of  $\text{CaCO}_3$  / L using Part 3's method.

#### ***Experimental***

Part 1 of the lab was dedicated to standardizing a solution of sulfuric acid. This was done via the process of titration. First, 1 gram of  $\text{Na}_2\text{CO}_3$  was measured and added to a 250mL flask and 50mL of distilled water to dissolve the salt. Another 50mL of distilled water was added to create the "Standard  $\text{Na}_2\text{CO}_3$  Solution," then its molarity was calculated using *Equation 1*.

*Equation 1: Used to find the molarity of Standard  $\text{Na}_2\text{CO}_3$  Solution*

$$\frac{1.03 \text{ grams of Na}_2\text{CO}_3}{0.200 \text{ L}} * \frac{99.96 \text{ grams of Na}_2\text{CO}_3}{100 \text{ grams}} * \frac{1 \text{ mol of Na}_2\text{CO}_3}{105.96 \text{ grams of Na}_2\text{CO}_3} = 0.0486 \text{ M of Na}_2\text{CO}_3$$

Next, the titration buret was conditioned by passing through sulfuric acid then was filled to the 0 mL line. Afterward, 20mL of the Standard  $\text{Na}_2\text{CO}_3$  Solution was added to a flask to begin the process of titration. To that flask, 20mL of distilled water and 3-4 drops of phenolphthalein were added. Before titration, this solution was a vibrant pink color. Next, the value for the first endpoint,  $EP_1$ , was calculated using *Equation 2*.

*Equation 2: Estimating the amount of sulfuric acid to reach  $EP_1$  (theoretical value)*

$$0.02 \text{ L} * \frac{0.0486 \text{ M}}{1 \text{ L of soln}} * \frac{1 \text{ mol of } H^+}{1 \text{ mol of } CO_3^{2-}} * \frac{1 \text{ mol of } H_2SO_4}{2 \text{ mols of } H^+} * \frac{1 \text{ L of soln}}{0.05 \text{ moles of } H_2SO_4} = 9.7 \text{ L of } H_2SO_4$$

Once the theoretical value was attained, the solution was titrated until it became a transparent color. The amount of sulfuric acid required to reach that clear color was recorded, then 3-4 drops of a mixed indicator were added to create a blue-green color. The flask was placed on a heat plate then cooled down in ice water. The solution was then titrated until the discernable blue-green color faded away. The amount of sulfuric acid required to achieve this transparency was noted as  $EP_2$ . The molarity of the sulfuric acid was calculated using *Equation 3*. This process has been repeated a total of three times to get three different trials. The results for this first part are recorded in *Results Table 1*.

*Equation 3: Used to find the molarity of sulfuric acid*

$$\frac{20 \text{ mL of standard soln}}{EP1 \text{ or } EP2} * \frac{0.0486 \text{ mol of } CO_3^{2-}}{1 \text{ L standard soln}} * \frac{1 \text{ mol of } H^+}{1 \text{ mol of } CO_3^{2-}} * \frac{1 \text{ mol of } H_2SO_4}{2 \text{ mols of } H^+} = M \text{ of } H_2SO_4 \text{ from } EP1 \text{ or } EP2$$

Part 2 of the experiment was to determine the alkalinity of an unknown. This experiment begins with the selection of one of many solutions: solution  $\epsilon$  for this experiment. 20 mL of the unknown solution, 20mL of distilled water, and 3-4 drops of phenolphthalein were added to a total of three flasks. Once added, the solution became a vibrant pink, which indicated a need for titration. The amount of mL required to titrate the solution was recorded as  $EP_1$ . Once titrated to

a clear solution, 3-4 drops of a mixed indicator were added. Once added, the solution became a blue-green color. This indicated that another titration process needed to occur. The flask was heated up, then cooled down and titrated until the blue-green color began to fade. The mL of sulfuric acid needed to titrate the solution was recorded as  $EP_2$ . Once every trial was completed, the alkalinities were calculated using *Equation 4*, with the  $A$  value being represented from *Results Table 3*.

*Equation 4: Used to calculate the alkalinities of the solution with given endpoint values*

$$\frac{A \text{ mL of } H_2SO_4}{20 \text{ mL sample soln}} * \frac{0.05 \text{ mol of } CaCO_3}{1 \text{ L of } H_2SO_4} * \frac{1 \text{ mol of } CaCO_3}{1 \text{ mol of } H_2SO_4} * \frac{100 \text{ g of } CaCO_3}{1 \text{ mol of } CaCO_3} * \frac{1 \text{ mg of } CaCO_3}{0.001 \text{ g of } CaCO_3} = \text{mg of } CaCO_3/L$$

Part 3 of this experiment was to estimate the potentiometric titration of an unknown. To begin, a pH meter was calibrated, then 20 mL of solution  $\epsilon$  and 50mL of distilled water were added to a beaker. The pH of the solution was measured in 0.5mL increments of the addition of  $H_2SO_4$ . The second set of points was calculated to determine the endpoints to determine the alkalinities using *Equation 3* and *Results Table 2*.

## Results and Discussion

Table 1: The average molarity of sulfuric acid and trials for Experiment 1

	<b>Trial 1</b>		<b>Trial 2</b>		<b>Trial 3</b>	
	$EP_1$	$EP_2$	$EP_1$	$EP_2$	$EP_1$	$EP_2$
Total vol of $Na_2CO_3$	20 mL	20 mL	20 mL	20 mL	20 mL	20 mL
Init vol of $H_2SO_4$	0 mL	10.6 mL	10 mL	20.8 mL	32 mL	35 mL
Final vol of $H_2SO_4$	10.6 mL	20.6 mL	20.8 mL	30.4 mL	42.2 mL	44.6 mL
Total vol of $H_2SO_4$	10.6 mL	10 mL	10.8 mL	9.6 mL	10.2 mL	9.6 mL
Molarity of $H_2SO_4$	0.0485 M	0.0486 M	0.0450 M	0.0506 M	0.0477 M	0.0506 M
Average M of $H_2SO_4$	0.0485 M					
Standard Deviation	0.0019					
% RSD	3.912 %					

Table 2: Calculations for A values to be used in Alkalinity

	Total Alk	Pheno Alk	Carb Alk	Bicarb Alk	Hydrox Alk
A Value	$EP_1 + EP_2$	$EP_1$	$2 * EP_1$	$EP_2 - EP_1$	0

Table 3: Titration data for solution  $\varepsilon$

	<b>Trial 1</b>		<b>Trial 2</b>		<b>Trial 3</b>	
	$EP_1$	$EP_2$	$EP_1$	$EP_2$	$EP_1$	$EP_2$
<i>Vol of the unknown sample</i>	20 mL	20 mL	20 mL	20 mL	20 mL	20 mL
<i>Init vol of <math>H_2SO_4</math></i>	0 mL	4.9 mL	9 mL	13.6 mL	12 mL	16.5 mL
<i>Final vol of <math>H_2SO_4</math></i>	4.9 mL	12.1 mL	13.6 mL	19.9 mL	16.5 mL	23 mL
<i>Total vol of <math>H_2SO_4</math></i>	4.9 mL	7.2 mL	4.6 mL	4.5 mL	4.5 mL	6.5 mL
<i>Total Alkalinity</i>	3025		2725		2725	
<i>Phenolphthalein Alkalinity</i>	1225		1150		1150	
<i>Carbonate Alkalinity</i>	2450		2300		2300	
<i>Bicarbonate Alkalinity</i>	575		425		425	
<i>Hydroxide Alkalinity</i>	0		0		0	

Table 4: Averages of Alkalinities from Table 2 for convenience

<i>Total Alkalinity</i>	2833.33 mg of $CaCO_3$ / L
<i>Phenolphthalein Alkalinity</i>	1166.67 mg of $CaCO_3$ / L
<i>Carbonate Alkalinity</i>	2333.33 mg of $CaCO_3$ / L
<i>Bicarbonate Alkalinity</i>	500 mg of $CaCO_3$ / L
<i>Hydroxide Alkalinity</i>	0 mg of $CaCO_3$ / L

Table 5: Plot points for pH of the solution

Line 1		Line 2	
$H_2SO_4$ mL	pH	$H_2SO_4$ mL avg	$-\Delta pH / \Delta mL$
0	9.81	0.25	0.36
0.5	9.63	0.75	0.38
1	9.44	1.25	0.34
1.5	9.27	1.75	0.5
2	9.02	2.25	0.84
2.5	8.6	2.75	1.78
3	7.71	3.25	1.36
3.5	7.03	3.75	0.52
4	6.77	4.25	0.44
4.5	6.55	4.75	0.4
5	6.35	5.25	0.34
5.5	6.18	5.75	0.38
6	5.99	6.25	0.22
6.5	5.88	6.75	0.28
7	5.74	7.25	0.38
7.5	5.55	7.75	0.3
8	5.4	8.25	0.42
8.5	5.19	8.75	0.76
9	4.81	9.25	1.86
9.5	3.88	9.75	2.22
10	2.77	10.25	1.92

Graph 1: The potentiometric titration of solution  $\epsilon$  made from graphing Table 5




Table 6: Average alkalinities and endpoint values from Graph 1

Endpoint 1	2.75 mL
Endpoint 2	9.75 mL
Total Alkalinity	3125 mg of $\text{CaCO}_3$ / L
Phenolphthalein Alkalinity	687.5 mg of $\text{CaCO}_3$ / L
Carbonate Alkalinity	1375 mg of $\text{CaCO}_3$ / L
Bicarbonate Alkalinity	1750 mg of $\text{CaCO}_3$ / L
Hydroxide Alkalinity	0 mg of $\text{CaCO}_3$ / L

The main differences between the data from *Table 4* and *Table 6* reside in the method used to attain the values. The data from Table 6 would be more accurate in terms of accuracy because the solution's pH was being taken in shorter intervals rather than being calculated from experimental data. Possible sources of error would include titrating too fast, pH probe not being calibrated correctly, and transfer errors.